

Ubuntu ROS Motion-Planning Computer

Neuro MoveIt! Motion Package

Robot State Publisher

- updateJointState
- updateCanulaLength
- updateProbeElementLength
- updateProbeTipLength

Planning Scene Manager

- updateScannerTransform
- sendScannerToScene
- removeScannerFromScene
- updateSkullMesh
- sendSkullToPlanner
- removeSkulFromPlanner
- updateEntryPointTransform
- updateTargetPointTransform
- displaySceneInRViz

Motion Planner

- computeWorkspace
- visualizeWorkspaceInRViz
- computeTrajectory
- transformTrajectorySpace
- previewTrajectoryInRViz

ROS-IGTLink Bridge Data Types

Bridge Params

- openIGTLink Server IP
- openIGTLink Server Port #

Incoming Topics (units: meters)

Transform /mri_scanner_to_robot

Polydata /skull_mesh

Point /entry_point or PolyData /entry

String /canula_length

String /probe_element_length

String /probe_tip_length

Point /target_point

Note: Eventually, this may be polydata of zone for ablation

String /service/compute_workspace

String /service/plan_to_target_point

String /service/plan_to_home_pose

Vector /robot_joint_state

Outgoing Topics (units: meters)

Polydata /workspace_poly

Vector /trajectories

Need to decide either entry point or entry zone

Slicer UI Computer (OS Agnostic)

Motion Planning Module for 3D Slicer

Head Segmenter

- addTissueImagesToVolume
- exportAsMesh

Robot Joint State

- sendJointStateUpdate

Workspace Visualization

- displayWorkspace
- specifyEntryPoint
- specifyCanulaOptions

Trajectory Control

- computeTrajectory
- showTrajectories
- chooseTrajectory
- executeTrajectory

Robot Controller

Controller (units: mm and degrees)

- executeTrajectory
- getJointState

Do we want to see trajectory live? We probably don't even need multiple options.

ROS

ROS

ROS

OpenIGTLink

OpenIGTLink